

CSG Dicoms Anonymizer



Stephen Larroque
Coma Science Group, GIGA research
University of Liège



07/04/2017



Advantages of CSG Dicoms Anonymizer

2

- ❑ **Automatically** detect name and **remove** it from any field, even **hidden ones**
- ❑ **Uniformize names** (invariant to typos & words switching)
- ❑ Can anonymize **demographics** along (same anon ids)
- ❑ Can use same demographics file for any set of dicom → anonymization will shorten to pertinent subjects
- ❑ Can continue a partial anonymization (eg, bug, access denied, ...)
- ❑ Deterministic anonymization → Can **update** anonymized demographics
- ❑ Generate list of **missing** demographics (ie, dicom files are present but no demographics for them).
- ❑ Generate a set of csv files to **deanonymize**.
Tip: These files can be encrypted in a 7z (not zip) archive with password and sent to the collaborator, he'll send it back if need more infos (ie, less work and files storage for us).

Dicoms Anonymizer – Algorithm

3

1. Generate list of patient names from dicoms (folders and zips)
2. Generate unique list of names (disambiguate similar names)
→ dicom_names.csv
3. Generate MD5 hash from names (with salt if provided → each lab can generate unique deterministic ids by tweaking the salt)
4. anonymized id = Shortened MD5 hash
→ idtonames.csv
5. If demographics: merge with dicoms names (compute distance matrix using disambiguation based on letters + words normalized levenshtein distance)
6. Apply anonymized id to dicoms files, folders names and demographics → anonymized dicoms & demographics
7. Delete non dicom files (pdf, doc, docx, txt, etc.)

Dicoms Anonymizer – Usage

4

1. Copy **all dicom** folders/zips in one folder
2. Get **demographics** file (optional)
3. Open **csg-fileutil dicoms anonymizer** (using Jupyter Notebook).
4. **Replace parameters** (dicom rootpath, demographics file path) **under each Part x**.
5. Click **Kernel > Restart & Run All**

→ Result: anonymized dicoms & demographics

Dicoms Anonymization – GUI

csg_fileutil_dicoms_anonymizer.exe

Settings
Dicoms anonymizer v1.3.2
Description: Anonymize dicoms and demographics using undecryptable hashes.
Note: use --cmd to avoid launching the graphical interface and use as a commandline tool.

Required Arguments


Input
Path to the dicom root folder (dicoms will be replaced and non-dicom files deleted).

Demographics
Demographics file to anonymize.

Optional Arguments

Distance Distance threshold for the jaccard distance (words <input type="text" value="0.2"/> <input type="button" value="Browse"/>	Anon_Prefix Prefix of the generated anonymized ids. <input type="text" value="subj_"/>
--	---

Dicoms Anonymization – Jupyter Notebook (old, please use GUI now)



localhost:8888/notebooks/csg_fileutil_dicoms_anonymizer.ipynb

jupyter csg_fileutil_dicoms_anonymizer Last Checkpoint: 03/18/2017 (autosaved)

File Edit View Insert Cell Kernel Widgets Help Python [default]

In []:

```
# Dicoms anonymizer
# By Stephen Larroque @ Coma Science Group, GIGA Research, University of Liege
# Creation date: 2017-02-07
# License: MIT
# v1.1.5
#
# INSTALL NOTE:
# You need to pip install pandas before launching this script.
# Tested on Python 2.7.11
#
# TODO:
# * Nothing critical.
#
```

In []:

```
# Forcefully autoreload all python modules
%load_ext autoreload
%autoreload 2
```

AUX FUNCS

In []:

```
# Auxiliary libraries and necessary functions

import os
import re
import shutil
import sys
```

Anonymized demographics

7

1	name	gender	age	final_dia	mri_sedation	accident_date	accident_etiology
2		M	54.0	MCS-	yes	25/04/2014	post arret cardiaque
3		M	57.0	EMCS	yes	25/08/2011	post traumatisme (c
4		M	49.0	UWS	yes	22/04/2003	post-anoxie (infarct
5		M	21.0	MCS+		26/10/2010	post trauma (le
6		M	19.0	MCS+	yes	30/07/2014	post trauma (le
7		F	46.0	UWS	no	5/01/2016	post-arret cardioresp
8		F	63.0	MCS+	yes	12/07/2010	post-avc ishemique



1	name	gender	age	final_dia	mri_seda	accident_	accident_	acquisiti
2	subj081	M	49.0	UWS	yes	#####	post-anox	15/03/201
3	subj084	F	46.0	UWS	no	#####	post-arret	31/05/201
4	subj086	F	51.0	coma	no	#####	post traum	02/12 - 17,
5	subj085	F	66.0	UWS	yes	#####	post hema	07/05/201
6	subj115	F	34.0	MCS+	no	#####	post-traur	18/04 - 23,
7	subj019	M	27.0		no			30/06/200
8	subj014	M	73.0	UWS	no	#####	accident :	12/07 - 19,

Resulting files

8

What you can **send**:

- Anonymized dicoms
- Anonymized demographics
(demographics_anonymized_shortened.csv)
- Missing demographics anonymized
(missing_demo_anonymized.csv)

What you need to **keep** (but not send):

- **idtoname.csv** → **anonymization mapping**, if collaborator might need more info about 1 subject
- dicom_names.csv → to regen anon (eg, to update demo)
- missing_demo.csv → missing demographics
- demographics_shortened.csv (optional)

Dicom Anonymization – Tips & tricks

9

- ❑ **Dicoms will be replaced**, advised to **backup (zip)** before anonymization (in case something went wrong and you need to restart)
- ❑ Anonymization can be continued if error or stopped (but disadvised)
- ❑ Demographics automatically **shortened** to subjects present in dicoms → can use the same demographics for all anonymizations
- ❑ Script divided in **3 independent parts**: 1. extract dicom names, 2. generate anonymization mapping, 3. anonymize dicoms & demographics.
→ **Can restart at any part (eg, to update demo)**

Thank you for your attention



James S. McDonnell Foundation



Université de Liège



BONUS SLIDES



Dicoms Anonymizer – Old Algorithm

12

1. Generate list of patient names from dicoms (folders and zips)
 2. Generate unique list of names (disambiguate similar names)
→ dicom_names.csv
 3. Generate MD5 hash from names
 4. Reorder names by MD5 hash
 5. New order = anonymized id
→ idtonames.csv
 6. If demographics: merge with dicoms names (compute distance matrix using disambiguation based on letters + words normalized levenshtein distance)
 7. Apply anonymized id to dicoms files, folders names and demographics → anonymized dicoms & demographics
 8. Delete non dicom files (pdf, doc, docx, txt, etc.)
- Future: add salt for anonymized id (so that each lab can generate its own unique deterministic ids)